


# Novel microfluidic developments for biomedical research applications



Alyssa C. Henry, Nicole Y. Morgan, Eugene J. Choi, Caitlyn M. McCullough, Seema Kacker, Elizabeth B. Brokaw, Paul D. Smith

## Microfluidic Devices: An Overview

We present a number of current developments in microfluidic systems for biomedical applications. The small length scale of microfluidic systems permits the use of smaller sample sizes, leading to decreased reagent costs, as well as faster analysis times. Several commercial microfluidic products have recently become available; however, there are ample opportunities in the biomedical research environment for improvements to existing technology as well as specialized applications.



Microfluidic devices can incorporate extraction, purification, amplification, labeling, separation, and detection of analytes in a complex matrix all on one small platform.

## Advantages of Microfluidics

- Sub-microliter sample volumes
- Rapid separation and analysis
- Numerous analytical steps can be performed in one small-footprint device
- Potentially less expensive than macro-sized counterparts:
  - less reagent consumption
  - multiplexing and mass-production easier

## Materials for Microfluidic Devices


Silicon/glass  
Chemically and physically robust  
Directly patterned with photolithography

Poly(dimethylsiloxane) (PDMS)  
Transparent elastomer  
Good for rapid prototyping  
Molded onto silicon templates

Polymers  
Hot imprinted or injection molded  
Facile thermal bonding  
Good for single-use devices

## Technological challenges


### Controlling fluid flow



**Electroosmosis** is the primary method of flow control in microchannels  
Requires charges on channel walls

### Developing Alternate Detection Schemes

Currently, most microfluidic systems use fluorescence for detection  
• requires fluorescent labeling of analytes  
• no information about molecular structure




## Photodirected Tethering of Charged Species in Microfluidic Devices for Fluid Flow Control


In making multifunctional devices, there is a need to create spatial variations of the surface properties along the length of a sealed microchannel. The ability to tailor these properties could lead to a variety of applications, such as improved electroosmotic pumping from highly charged channel regions, or the targeted patterning of capture ligands. In the present study, we wish to tether moieties with an additional charged species to the microchannel.

### Detailed chemistry


1. React microchannel surface with a silane to yield a reactive functional group.




2. Attach photoprotecting group: 4,5-dimethoxy-2-nitrobenzyl carbonyl (NVOC)




4. Attach carboxylated poly(ethylene glycol) (PEG) to minimize non-specific adsorption and to provide strongly negative surfaces for EOF



### Patterning inside a microchannel: Monolayer photolithography




3. Selectively expose to 355 nm light to get spatially patterned reactive primary amines.




5. After tethering PEG to sections of the microchannel, deprotect other portions of the microchannel for further reactions.


### Contact angle goniometry




### Fluorescence microscopy studies of glass surfaces terminated in various organic moieties



### UV-visible spectroscopy studies of glass surfaces terminated in various organic moieties



### Electroosmotic flow studies




## Surface-Enhance Raman Spectroscopy as a Detection Method in Microfluidic Devices


Currently, the most common detection method used in microfluidic systems is fluorescence detection. While fluorescence offers high sensitivity, chemical labeling of analytes with fluorophores is generally required. In addition, non-specific labeling and autofluorescence is often a problem. We present recent efforts to use patterned Au colloids in a microdevice for surface-enhanced Raman spectroscopy (SERS) detection. This detection technique offers direct measurement of biomolecules in aqueous solution.

### Synthesis and characterization of Au colloids and Au colloid-terminated glass


Au colloids were synthesized using  $\text{HAuCl}_4$  and sodium citrate. The resulting citrate-stabilized Au colloids were 13 nm in diameter, with a negatively-charged shell.



### UV-Visible absorption of Au colloids in solution




### UV-Visible absorption spectra of Au colloids adsorbed on amine-terminated glass surfaces




- Au colloids in solution have an absorption maximum at ~520 nm
- Au colloids adsorbed on amine-terminated surfaces have an absorption maximum similar to the colloids in solution
- The absorption maximum of the Au colloid monolayers shifts after immersion in the Au colloids solution for long periods

### Synthesis of Au colloid multilayers




### Schematic of Optical Detection System



## Rapid tuning of protein deposition conditions for AFM

Elastomeric microchannels are used to temporarily create a microdevice on an atomic force microscopy (AFM) substrate. Because flow in microchannels is laminar, all mixing is diffusive, which allows us to create a well-controlled buffer gradient, either in ionic strength or pH, over a few hundred microns. This permits a rapid sampling of buffer conditions for protein deposition on AFM substrates.




Top view: two side arms used as inlets for protein solutions in different buffers

First, the top of cross is used as the outlet, in order to balance flow from the sides before exposing lower part of device to protein.

Target area for imaging is in bottom of cross




### Device Operation: Salt concentration gradient




### Cytochrome C: salt gradient

As the concentration of NaCl increases across the channel, the protein deposition decreases.



### Device Operation: pH gradient



Bovine Serum Albumin: pH gradient  
66kDa, ~4x8nm  
Isoelectric point: 4.9

